

Identification of a Novel *Moraxella* Species in a Culture-Negative CNS Shunt Infection Using Metagenomic Next-Generation Sequencing

Denver Niles^{1,2}, Sirisha Thippabhotla³, Nick Gimbrone³, Brian O'Donovan³, Sudeb Dalai³, Steve Miller³

1. Baylor College of Medicine
2. Texas Children's Hospital, Houston, TX, USA
3. Delve Bio, San Francisco, CA, USA




Background

- Metagenomic next-generation sequencing (mNGS) enables rapid, comprehensive, and culture-independent detection of pathogens, including atypical and novel species, directly from CSF, and has demonstrated clinical utility in challenging clinical scenarios.
- CNS shunt infections are often culture-negative due to fastidious pathogens, low microbial burden, or prior antimicrobial treatment.
- A case of a pediatric CNS shunt infection is described in which CSF gram stain was positive but multiple cultures and other tests were repeatedly negative; mNGS ultimately identified a potentially novel *Moraxella* species.

mNGS Analysis


Automated protocol focused on precision, sterility, speed, and throughput



- Test: Delve Detect CSF (Delve Bio), CLIA-certified/CAP-accredited laboratory- developed test (LDT)
- Sample: 1 mL CSF, DNA/RNA deep sequenced (>10M reads/library)
- Findings:

- Moraxella* bacteria was detected at subthreshold levels, despite the presence of high host background
- Sequences aligned to multiple *Moraxella* species, making species-level identification challenging
- Genomic coverage over the ribosomal RNA operon (RRN) was sufficient to build a consensus sequence and perform phylogenetic analysis
- The closest phylogenetic matches by RRN were *Moraxella* sp. K1664 along with other Group I *Moraxella* species
- Phylogenetic placement between these multiple species indicates a potentially novel *Moraxella* organism

Findings



Subthreshold Positive Result

- Moraxella* species detected (below reporting threshold), closely aligning with *Moraxella bovis* and atypical species.
- Sample had high host DNA background, reducing sensitivity for detecting pathogens.




Figure 1. Elevated *Moraxella* Biomass in Patient Sample Compared to Historical mNGS Samples

Histogram of \log_{10} -transformed total microbial biomass (pg) for *Moraxella* species detected in CSF samples submitted for clinical mNGS testing. The patient sample (red arrow) demonstrates significantly higher organism biomass than the background level seen in historical clinical samples, supporting the presence of true infection rather than environmental contamination or incidental signal.




Figure 2. Read Alignment to *Moraxella* sp. K1664 Genome

Genome coverage plot showing mNGS read alignment to *Moraxella* sp. K1664 reference genome. Upper panel shows genome-wide consensus identity and read coverage at 1.7% of the genome. Lower panel shows zoomed-in coverage across the ribosomal RNA operon (RRN) which was used for phylogenetic analysis.

Species	rpM	rpM Ratio	Unique Reads	Biomass (fg)
<i>Moraxella</i> (declassified reads)	0.22991	0.22991	8	42.123
<i>Moraxella bovis</i>	2.21291	2.21291	66	405.434
<i>Moraxella nonliquefaciens</i>	0.31613	0.31613	9	57.9191
<i>Moraxella bovovuli</i>	0.17244	0.17244	6	31.5922
<i>Moraxella haemolytica</i>	0.08622	0.08622	1	15.7961
<i>Moraxella ovis</i>	0.02874	0.02874	1	5.26537
<i>Moraxella nasibovis</i>	0.02874	0.02874	1	5.26537
<i>Moraxella</i> genus (total)	3.07509	3.07509	92	563.39514

Table 1. Sequencing Metrics for *Moraxella* DNA Library Reads Detected by mNGS

Table shows reads per million (rpM), rpM ratio (normalized abundance relative to non-template control signal), number of unique reads (by exact sequence match), and microbial biomass (fg) for individual *Moraxella* species observed. *Moraxella* (species-declassified reads) indicates reads assigned to the genus level without confident species-level classification. *Moraxella* genus (total) indicates the sum across all *Moraxella* species and declassified reads. Elevated rpM and biomass were observed for *M. bovis* with lower contributions from several other *Moraxella* species, supporting the presence of a novel or divergent organism within the *Moraxella* genus.




Figure 3. Phylogenetic Analysis of *Moraxella* Consensus Sequence

Maximum likelihood phylogenetic tree constructed from patient-derived ribosomal RNA operon (RRN) consensus sequence demonstrates that the patient sample (blue) clusters within Group I *Moraxella*. The patient sequence appears closely related to *Moraxella* sp. K1664, with approximately equivalent distance between several closely-related *Moraxella* species, supporting the presence of a potentially novel *Moraxella* species. (Scale: Tamura-Nei distance = 0.01)

Discussion

- Despite antibiotic pre-treatment, mNGS identified a pathogen consistent with CSF Gram stain result in a case repeatedly negative by culture and multiple other tests
- mNGS provided definitive diagnosis and could have enabled narrowing of antibiotics
- CNS shunt-associated infections often involve fastidious or biofilm-forming organisms that can evade culture, particularly with prior antimicrobial exposure
- Moraxella* is an uncommon CNS pathogen and rarely reported in shunt infections; its identification here underscores:
 - Diagnostic limits of culture-based microbiology
 - mNGS value in high-risk, culture-negative infections
- Phylogenetic analysis supports a potentially uncharacterized *Moraxella* species.

Conclusions

- Delve Detect mNGS enabled the diagnosis of a culture-negative CNS shunt infection caused by a potentially novel *Moraxella* species.
- This case underscores the clinical utility of mNGS in complex infections and the ability to detect emerging or atypical pathogens missed by conventional diagnostics.
- Together with findings from other studies, this case supports broader integration of mNGS into the diagnostic workup for unexplained or culture-negative CNS infections, particularly in scenarios of CSF-associated shunts, devices, or other hardware